


HOSPITAL  
INNOVATIONS  
Would like to  
Introduce

## PULLUP® & PULLUP® XL **Stability** and **Choice**



| FEATURE                                      | PULLUP® | FIXED LOOP | ADJUSTABLE LOOP |
|--|---------|------------|-----------------|
| Proven fixation strength                     | ✓       | ✓          | ✗               |
| No clinically significant graft displacement | ✓       | ✓          | ✗               |
| Intra-operative flexibility                  | ✓       | ✗          | ✓               |
| Reliable flip                                | ✓       | ✓          | ✗               |

## No **Compromises**

**Next generation** adjustable loop suspensory fixation

## ■ ■ Confident Stability

Novel, patented locking mechanism designed to minimise graft slippage during tensioning<sup>1</sup>

Secondary positive flip suture ensures accurate seating of the button against the cortex<sup>1</sup>

## ■ ■ Intra-operative Choice

PULLUP® adjustable loop eliminates the need for intra-operative calculations and avoids limitations in tunnel and graft length choice<sup>1</sup>

PULLUP® fills the tunnel with tendon graft, avoiding concerns on graft stability and bone-tendon healing<sup>1</sup>

### Unique Loop-Locking System

Chinese finger trap configuration with simultaneous tensioning of loop sutures minimises loop lengthening and graft displacement<sup>1</sup>

### Adjustable Loop System

Allows for surgeon flexibility and eliminates the need to plan intra-operatively<sup>1</sup>

### Positive Suture Flip

Separate flip suture ensures reliable flipping, placing the button securely against the cortex<sup>1</sup>